

Treatment of the enuresis (bedwetting)

Technological innovation

1932 – 2017
85 years Pipi-Stop®

Vibrating wristband VIBRA wireless...

...an accessory to the Pipi-Stop® (basic device)

Discover our last development:

The vibrating wristband to use with or without the acoustic alarm.



Swiss made

The vibrating wristband VIBRA is worn on the wrist or on the arm

Most of the children who wet their bed sleep very deeply and are not aware of their nocturnal urination. Therefore it is **wise to begin the treatment with our Pipi-Stop which has an acoustic alarm**. This is the **best way to inform the parents** about the urination moment, so that they can wake the child up. In this way, the child will develop the reflex to retain his pee at the right moment.

When the child wakes up easier during the night, you can do without the acoustic alarm and continue with the vibrating wristband. This will avoid waking up other members of the family.


The **vibrating wristband** is used as an **accessory to the Pipi-Stop (basic device)**.

1st phase of the treatment:

At the beginning of the treatment it is necessary to use the Pipi-Stop (basic device) with the acoustic alarm.



The detector (conductive sheet of polyester) is placed in the underwear. In this way the detector stays in place. At the first drops of urine the transmitter sends a signal to the receiver: the alarm is activated (sound + flashlight).



2nd phase of the treatment:

As soon as the child starts waking up by himself at the sound of the acoustic alarm, it is possible to switch to vibrating mode (the wristband VIBRA can be used **with or without the acoustic alarm** of the receiver).

