

# Bedwetting

The nocturnal enuresis

## Pipi-Stop®

electronic alarm-system wireless  
rented for enuretic children



Swiss made

## Bedwetting!

This problem is much more frequent than we imagine. Even so, the subject is **taboo**. Not only the child who wets his bed (enuretic child) suffers hardly but also the whole family gets under pressure through this.

About 10% of the 6 years old children and about 6% of the 10 years old children wet their bed while they sleep. Different studies show that the enuresis is hereditary for 70 to 85% of enuretic children who are more than 5 years old, and boys are more frequently affected than girls.

These figures prove that **you are not the only one concerned by this problem**. There is no reason to over-dramatize. **Your child can also be cured from this illness and regain self-confidence.**

Generally, children stop wetting their bed before the age of 5 or 6 years. Some may need more time. It depends on the maturity of the child and of his physical development.

## Has the bedwetting become a handicap for your child?

Can he not take part any more to certain activities like other children (sleep at a friend's place, stay in a school camp, etc)?

If it is so, **it is the right moment to take this problem very seriously!** Together we will do everything to eliminate this inconvenience.

Our method is based on the **conditioning of the Pavlov reflex** (behavioural therapy of the enuresis nocturna).

**Pipi-Stop** is a healing method **without**

**medicine**, in the form of an **electronic alarm-system**.


Most of the children who wet their bed sleep very deeply and therefore are not aware of their nocturnal involuntary urination.

Through the alarm which will always be activated at the right moment (at the first drops of urine), the child will be able to learn to react and to **develop subconsciously the conditioned reflex to retain his pee**. If his bladder is full, he will be able to wake up and to go to toilet.

## Pipi-Stop is very simple to use

The detector (conductive sheet of polyester) is placed in the underwear. In this way the detector stays in place.

At the first drops of urine the transmitter sends a signal to the receiver: the alarm is activated (sound + flashlight).



## For a successful treatment

- The child must be motivated and participate to his therapy
- During the alarm the child must be **completely woken up** (conscious). It is necessary that an adult helps him to stand up, speaks to him, so that he remembers that the next morning
- The **child should be woken up only once in the night** (the learning effect is the same)
- Before you begin the treatment, it is better to consult a doctor to exclude any anatomic/biological cause of the enuresis.

## The advantages with Pipi-Stop

- ✓ The **most effective way** for the treatment of the enuresis (~ 90 % of success)
- ✓ The way that **counts the least of relapses**
- ✓ **No chemical substances**, therefore **no secondary effects**
- ✓ **Wireless**: the child is free in his movements. No risks to get tangled!
- ✓ **No danger**: the emitted waves are very weak (more than 100'000 times weaker than those of a mobile phone) and are only transmitted while the alarm rings
- ✓ **No drinking restriction**, even at the evening
- ✓ **Relatively short treatment**: the average is between 6 and 12 weeks
- ✓ **Very long experience**: since more than **85 years** we treat successfully the children who wet their bed. Hundreds of testimonies and letters of gratitude confirm continuously the effectiveness of our method
- ✓ **Follow-up of the cure**: consultancy and information about the treatment are included in the renting price
- ✓ **Pipi-Stop** fulfils the RoHS standards (**CE standards**).

In order to profit from the **refunding by your basic insurance** the child must be **5 years of age**. Do not forget to ask a prescription from your GP/paediatrician.

**We rent out our Pipi-Stop electronic alarm-system.**

To order it or for any further information, please contact us or visit our website [www.pipi-stop.ch](http://www.pipi-stop.ch)

## Melebi SA

Parkweg 12  
2502 Biel/Bienne  
tel. 032 751 27 95  
fax 032 751 27 82

mail: [info@melebi.ch](mailto:info@melebi.ch)  
[www.pipi-stop.ch](http://www.pipi-stop.ch)



Swiss made

CE Our Pipi-Stop fulfils the RoHS standards