

Pipi-Stop Instruction Manual

Components of the Pipi-Stop® warning device



Urine Detector



Flannel Pouch



Receiver




Transmitter



Power Supply

TESTING THE SYSTEM

To view test videos: www.pipi-stop.ch/en/pipi-stop-therapy/instruction-for-use/



Testing the receiver:

1. Connect the power supply to a power outlet (230V) and connect to the back of the receiver
2. Briefly press **one of the 2 volume control buttons**
3. **Test result point 1 and 2:**
blue light on + sound

Testing the transmitter:

1. Bridge (connect) the 2 metal studs with a small metal spoon
2. **Test result:**
flashing white light
+ sound





Testing the complete system:

1. Attach the transmitter to the detector (*the detector must be unpacked*)
2. Bridge the tongue with a spoon
3. **Test result:** flashing white light + sound



Important!

Testing with water will not set off the alarm!

If the alarm does not go off: make sure that the receiver and transmitter are in close proximity. If the Pipi-Stop does not work after these checks, **please contact us immediately.**

Safety / Maintenance

The receiver and power supply must NEVER come into contact with water (danger of electric shock!)

The transmitter must NEVER be opened. The battery can only be changed by Melebi SA for watertightness purposes.

Use **ONLY** our 12V power supply with our receiver to avoid overvoltage!

Any necessary repairs are carried out **exclusively** by Melebi SA.

Preparing the transmitter and the urine detector

Attach the transmitter to the urine detector
(the polarity direction
does not matter)



①

②

③




Slide the urine detector into the flannel pouch so that the transmitter is facing upwards. Then fold the flap of the pouch towards the inside and ensure that everything is well fixed.

④

⑤

⑥



Then place the prepared pouch in a slip and make sure that the transmitter is on the outside (*not against the body*), with the skinny part between the legs.

Preparing the receiver

Connect the power supply to the back of the receiver and connect to a 230-volt outlet
(blue indicator light comes on + beep tone)

Place the receiver on the child's bedside table




Adjusting the volume

Increase the volume: press the button (↑) and hold until you reach the desired loudness

Decrease the volume: press the button (↓) and hold until you reach the desired loudness

To activate the **silent mode**, press both volume buttons at the same time.
To deactivate, simply press one of the two buttons.



After nighttime urination

To stop the alarm:

1. Disconnect the transmitter from the urine detector

*If the alarm was longer than **3 minutes**, press briefly one of the two volume buttons on the receiver to stop the alarm*

2. Do not disconnect the power cable
to stop the alarm !

Cleaning the urine detector, transmitter and flannel pouch:

Separate the individual parts and rinse thoroughly under lukewarm running water immediately.



Rinse the transmitter and urine detector and dry them carefully
(**DO NOT RUB THE DETECTOR!**)

After cleaning the flannel pouch,
wring it out, hang it up and let it air
dry. (quick-drying)

Melebi SA

Parkweg 12
2502 Biel/Bienne
phone 032 751 27 95
fax 032 751 27 82

email: info@melebi.ch
www.pipi-stop.ch



CE Our Pipi-Stop fulfills
RoHS standards